


Technical Drawing

Sealing Mats for 0.5 and 1.2 ml Deep Well Plates

Order number 96.09711.9.01




- Material: TPE
- Pierceable, Non-Autoclavable Mats
- Solvent resistant (including DMSO)
- Temperature range -20°C to + 55°C

Dimensions : Millimeters

Technical Drawing

Sealing Mats for 2.2 ml Deep Well Plates

Order number 96.09537.9.01



- Material: EVA
- Non-Autoclavable Mats
- Solvent resistant (including DMSO)
- Temperature range -20°C to + 55°C

Dimensions : Millimeters