



Single Channel Transferpette S®											Maximum Volume					Multi Channel Transferpette S®					Pipette Tips
1 µl	2.5 µl	10 µl	20 µl	25 µl	50 µl	100 µl	200 µl	1'000 µl	5 ml	10 ml	10 µl	50 µl	100 µl	200 µl	300 µl						
•	• <sup>2)</sup>	• <sup>2)</sup>									•					0.1 - 20 µl					
	•	•									•					0.5 - 20 µl					
	• <sup>2)</sup>	•									•					1 - 50 µl					
			•	•	•	•	•					•	•	•	• <sup>4)</sup>	2 - 200 µl <sup>1)</sup>					
			•	•	•	•	•					•	•	•	•	5 - 300 µl					
								•								50 - 1'000 µl <sup>1)</sup>					
									•							0.5 - 5 ml <sup>3)</sup>					
										•						1 - 10 ml <sup>3)</sup>					
•																0.1 - 1 µl					
	•	•									•					0.5 - 10 µl					
	•	•									•					1 - 20 µl					
			•													2 - 20 µl					
			•	•	•	•						•	•	• <sup>4)</sup>	• <sup>4)</sup>	5 - 100 µl					
			•				•					•	•	•	• <sup>4)</sup>	5 - 200 µl					
								•								50 - 1'000 µl					

- 1) Tips in bags are yellow or blue. Racked or stacked tips are clear and the carrier plate is yellow or blue.
- 2) Calibration or adjustment may be necessary
- 3) Filter for these tips are available separately
- 4) Max. volume of tip is below max. volume of pipette

