

Nolato Treff Quality

CleanRoom Pure[®]

Basic Quality and Purity Standard of all TreffLab Products

“Absolutely PURE Right From The Start”



Statement

Nolato Treff AG confirms that TreffLab consumables and products marked with the CleanRoom[®] Pure quality symbol have been manufactured in microbiologically controlled clean rooms class 8, ISO 14644-1. To avoid direct human contact, employees are fully protective dressed and masked.

Quality Parameters

Quality standard for original sealed packages


- Free of RNase (< 5x10E-9 Kunitz units)
- Free of DNase (< 5x10E-6 Kunitz units)
- Free of DNA (< 2 pg DNA)
- Free of endotoxins (< 10 copies)
- Free of PCR inhibitors (< 0.001 EU/ml)
- Free of R and ATP (< 30 RLU)

Analyses are performed by an independent Swiss laboratory for molecular biological & microbiological analytics, accredited ISO/IEC17025 for real-time PCR.

03. March, 2021



Richard Wiesli
Sales & Marketing Manager,
TreffLab



Thomas Kreil
Quality Department