


Nolato Treff Quality

LoTrace[®] DNA

**Low binding DNA
PCR Grade Surface**

**“Never lose your valuable
DNA trace sample”**



Statement

Nolato Treff AG confirms that TreffLab consumables and products marked with the LoTrace[®] DNA quality symbol have been manufactured in a special additive-free process to achieve Low Binding DNA – PCR Grade surface. All products are manufactured in microbiologically controlled clean rooms class 8, ISO 14644-1. To avoid direct human contact, employees are fully protective dressed and masked.

Quality Parameters

Quality standard for original sealed packages

- | | |
|--------------------------|--------------------------|
| ▶ Low Binding surface | Real-time PCR |
| ▶ Free of RNase | (< 5x10E-9 Kunitz units) |
| ▶ Free of DNase | (< 5x10E-6 Kunitz units) |
| ▶ Free of DNA | (< 2 pg DNA) |
| ▶ Free of endotoxins | (< 10 copies) |
| ▶ Free of PCR inhibitors | (< 0.001 EU/ml) |
| ▶ Free of R and ATP | (< 30 RLU) |

Analyses are performed by an independent Swiss laboratory for molecular biological & microbiological analytics, accredited ISO/IEC17025 for real-time PCR.

03. March, 2021

Richard Wiesli
Sales & Marketing Manager,
TreffLab

Thomas Kreil
Quality Department