



## Materials

All products are produced with 100 % virgin polymer only. No recycled materials are used. Our material suppliers do not use or intentionally incorporate the following agents into the materials

Nolato Treff used for production:

- Slip agents (i.e. Oleamide)
- Biocides (i.e. DIHEMDA)
- Plasticizers (softeners/ phthalates)
- Silicone

## ISO Certificates

The Nolato Treff AG is certified according to ISO 9001, ISO 13485, ISO 14001 and ISO 22000.

All TreffLab Consumables are produced in microbiologically controlled clean rooms class 8, ISO 14644-1

## Safety note

Additional sterilized Pipettes Tips are available on special request. Sterile pipette versions are sterilized with  $\beta$ -radiation according to ISO 11137 and AAMI guidelines. The sterility meets USP 29 and Ph. Eur. requirements SAL of  $10^{-6}$  ( $\beta$ -radiation). Intended use: For Research and Development only.


TreffLab Slim-Line™ Pipette Tips, PP, are designed to achieve the highest level of precision, accuracy and reproducibility demanded by scientists for challenging assays and analysis. To keep the physical properties of TreffLab PP Consumables, the following storage conditions are recommended:

- storage temperature (18 – 28°C)
- no light exposure (= packed in original cardboard box)

March, 2022



Richard Wiesli  
Sales & Marketing Manager,  
TreffLab



Thomas Kreil  
Quality Department

